



User's Manual

Symmetric Dimethylarginine ELISA Kit



DEIA093J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The SDMA ELISA Kit is intended for the quantitative determination of symmetric dimethylarginine (SDMA) in human EDTA-plasma and serum.

Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

Precision

Intra-assay Precision: 4.8%-7.5%

Inter-assay Precision: 6%-7%

Sensitivity

0.05 µmol/L