



User's Manual

Thyroglobulin recovery ELISA Kit



DEIA096J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The described ELISA is intended for the quantitative determination of human thyroglobulin (hTG) in plasma and serum.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Sensitivity

1 µg/L

Specificity

The described hTG assay is highly specific. No cross reactivity with other proteins in plasma and serum was observed.

