



User's Manual

Human Chymotrypsin ELISA Kit



DEIA10041



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Chymotrypsin ELISA Kit is a sandwich Enzyme Immuno Assay intended for the quantitative determination of Chymotrypsin in stool.

Storage

2-8°C

Precision

Intra-assay CV n= 20

Sample	Chymotrypsin, mean value [ng/ml]	Intra-Assay CV [%]
1	24,20	3,7
2	16,35	6,9

Table 3.

Inter-assay CV n= 31

Sample	Chymotrypsin, mean value [ng/ml]	Inter-Assay CV [%]
1	101,8	5,3

Table 4.

