



## User's Manual

# Rat Alarin ELISA Kit



**DEIA10446**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Rat Alarin ELISA Kit is designed to detect a specific peptide and its related peptides in rat Blood, Plasma, Tissue and CSF, based on the principle of "competitive" enzyme immunoassay.

### Storage

Store the kit at 4°C upon receipt.

### Detection Range

0-100 ng/mL

### Sensitivity

0.11 ng/mL