



User's Manual

Human, Porcine, Ovine, Canine ANP, alpha (1-28) ELISA Kit



DEIA10463



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human, Porcine, Ovine, Canine ANP, alpha (1-28) ELISA Kit is designed to detect a specific peptide and its related peptides in Human, Porcine, Ovine and Canine Blood, Plasma, Tissue and CSF, based on the principle of "competitive" enzyme immunoassay.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-25 ng/mL

Sensitivity

0.12 ng/mL