



**User's Manual**

# Human Free Beta-Subunit of Chorionic Gonadotropin ELISA Kit

REF

DEIA10588



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  **Web:** [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human Free Beta-Subunit of Chorionic Gonadotropin ELISA Kit is for the quantitative determination of free beta subunit of human chorionic gonadotropin (free Beta-hCG) concentration in human serum.

### Storage

Unopened test kits should be stored at 2-8°C upon receipt. For more detailed information, please download the following document on our website.

### Sensitivity

0.25 ng/mL