



**User's Manual**

# Human INSL5 C-Peptide, Prepro (49-106) ELISA Kit



DEIA10637



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Human INSL5 C-Peptide, Prepro (49-106) ELISA Kit is intended for the detection of insulin-like 5 in Blood, Plasma, Tissue, CSF.

### Storage

Store the kit at 4°C upon receipt.

### Detection Range

0.1-100 ng/mL

### Sensitivity

0.33 ng/mL