



User's Manual

Human INSL-7, C Peptide ELISA Kit



DEIA10640



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human INSL-7, C Peptide ELISA Kit is intended for the detection of relaxin 3 in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL