



User's Manual

Human LANP/ANP (26-55), prepro ELISA Kit



DEIA10649



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human LANP/ANP (26-55), prepro ELISA Kit is intended for the detection of acidic (leucine-rich) nuclear phosphoprotein 32 family, member A in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.32 ng/mL