



User's Manual

Human, Rat, Mouse Porcine LH- RH/Gn-RH ELISA Kit



DEIA10654



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse, Porcine LH-RH/Gn-RH ELISA Kit is intended for the detection of gonadotropin-releasing hormone 1 (luteinizing-releasing hormone) in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-25 ng/mL

Sensitivity

0.08 ng/mL