



User's Manual

Human, Rat, Mouse Neurokinin B/Neuromedin K ELISA Kit



DEIA10677



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse Neurokinin B/Neuromedin K ELISA Kit is designed to detect a specific peptide and its related peptides in Human, Rat and Mouse Blood, Plasma, Tissue and CSF, based on the principle of "competitive" enzyme immunoassay.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.14 ng/mL