



User's Manual

Human, Rat, Mouse Neuropeptide EI ELISA Kit



DEIA10683



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse Neuropeptide EI ELISA Kit is designed to detect a specific peptide and its related peptides in Human, Rat and Mouse Blood, Plasma, Tissue and CSF, based on the principle of "competitive" enzyme immunoassay.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.25 ng/mL

