



User's Manual

Human, Rat, Mouse, Porcine, Ovine, Bovine, Monkey Orexin A ELISA Kit



DEIA10700



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse, Porcine, Ovine, Bovine, Monkey Orexin A ELISA Kit is intended for the detection of Prepro-orexin in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.22 ng/mL