



User's Manual

Human, Rat, Mouse Oxyntomodulin ELISA Kit



DEIA10702



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse Oxyntomodulin ELISA Kit is intended for the detection of Oxyntomodulin in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.25 ng/mL