



User's Manual

Human Salusin-alpha ELISA Kit



DEIA10741



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human Salusin-alpha ELISA Kit is intended for the detection of Salusin-alpha in Blood, Plasma, Tissue, CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-100 ng/mL

Sensitivity

0.24 ng/mL