



User's Manual

Human, Rat, Mouse, Porcine, Ovine VIP ELISA Kit



DEIA10793



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human, Rat, Mouse, Porcine, Ovine VIP ELISA Kit is designed to detect a specific peptide and its related peptides in blood, plasma, tissue and CSF.

Storage

Store the kit at 4°C upon receipt.

Detection Range

0-25 ng/mL

Sensitivity

0.12 ng/mL