



User's Manual

Mouse C-MET/HGFR ELISA Kit



DEIA10993



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Mouse C-MET/HGFR ELISA Kit is for quantitative detection of mouse C-MET in cell culture supernates, serum and plasma (heparin, EDTA).

Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)

Detection Range

62.5 pg/mL-4000 pg/mL