



**User's Manual**

# **Spirometra mansonii IgG ELISA Kit**



**DEIA1102**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The test take the *Spirometra mansoni* eggs antigen as the solid phase antigen to detect IgG class antibodies against *Spirometra mansoni* in serum.

### Storage

Store at 4°C for frequent use, at -20°C for infrequent use. Avoid multiple freeze-thaw cycles. For more detailed information, please download the following document on our website.



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)