



## User's Manual

# Human Troponin T ELISA Kit



DEIA1102L



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human Troponin T ELISA kit is an in vitro enzyme-linked immunosorbent assay for the quantitative measurement of human Troponin T in serum (human Troponin T concentration is pretty low in normal serum/plasma and may not be detected in this assay), plasma, cell culture supernatants and urine.

### Storage

May be stored for up to 6 months at 2 to 8°C from the date of shipment. Opened Microplate Wells or reagents may be store for up to 1 month at 2 to 8°C. Return unused wells to the pouch containing desiccant pack, reseal along entire edge. Reconstituted standard can be stored at -80°C for up to 1 week.

### Sensitivity

0.35 ng/mL