



User's Manual

Human Procollagen I C- Terminal Propeptide(PICP) ELISA Kit

REF

DEIA1103L



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Kit is to be used for the vitro quantitative determination of Target protein in serum, plasma, tissue lysates, cell culture supernatants and other biological fluids. The Kit is intended for research use only, not for diagnostic or therapeutic procedure.

Storage

All the components in the kit should be stored up to 1 year at -20°C and 3 months at 2-8°C, and should be kept according to the labels on vials. If the kit can not be used within 2 weeks, please keep in -20°C.