



User's Manual

Human β Amyloid 42 antibody ELISA Kit



DEIA1128



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human β Amyloid 42 antibody ELISA kit is designed to detect and quantify the level of human β Amyloid 42 antibody in urine.

Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.