



User's Manual

Human Interleukin 12, IL12 ELISA Kit



DEIA117



10 plates



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Human IL-12 ELISA development kit contains the key components required for the quantitative measurement of natural and/or recombinant IL-12 in a sandwich ELISA format within the range of 32-3,000 pg/mL. Using the ELISA protocol described below, this kit provides sufficient reagents to assay IL-12 in approximately 1,000 ELISA plate wells.

Storage

Store at 2-8°C. For more detailed information, please download the following document on our website.



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com