



## User's Manual

# Human SSA Antibody ELISA Kit

REF

DEIA1220



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human SSA Antibody ELISA kit is designed to detect and quantify the level of human SSA Antibody in cell culture supernatant, serum, plasma and tissue.

### Storage

Store all reagents at 2-8°C. For more detailed information, please download the following document on our website.

### Sensitivity

65.0%