



**User's Manual**

# Human M2 type mitochondria Antibody ELISA Kit

REF

DEIA1224



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Human M2 type mitochondria Antibody ELISA kit is designed to detect and quantify the level of human M2 type mitochondria Antibody in cell culture supernatant, serum, plasma and tissue.

### Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

### Sensitivity

The sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is  $\geq 90\%$ .