



User's Manual

Human NUP210 Antibody ELISA Kit



DEIA1226



96T




This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human NUP210 Antibody ELISA kit is designed to detect and quantify the level of human NUP210 Antibody in cell culture supernatant, serum, plasma and tissue.

Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

The sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 29.0%.