



## User's Manual

# Human IL-12 p70 ELISA Kit



DEIA1335



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human IL-12 p70 ELISA kit is designed to detect and quantify the level of Human IL-12 p70 in cell culture supernatant, serum, and plasma.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Sensitivity

< 2.0 pg/mL