



User's Manual

Mouse IL-2 ELISA Kit



DEIA1349



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Mouse IL-2 ELISA kit is designed to detect and quantify the level of Mouse IL-2 in cell culture supernatant, serum, and plasma.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Sensitivity

< 1 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com