



User's Manual

Mouse IL-17E ELISA Kit

REF DEIA1354

 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Mouse KC(IL -8) ELISA kit is designed to detect and quantify the level of Mouse KC(IL -8) in cell culture supernatant, serum, and plasma.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.