



User's Manual

IL1RA antagonist Human ELISA Kit



DEIA1368



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human IL1RA antagonist concentrations in Serum, Plasma, Cell Culture Supernatants and Urine.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0.124 ng/mL-30 ng/mL

Sensitivity

0.1 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com