



User's Manual

IL2 Human ELISA Kit



DEIA1372



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human IL2 ELISA is to be used for the quantitative determination of interleukin-2 (IL2) in Human serum, plasma, buffered solutions or cell culture medium. The assay will recognize both natural and recombinant Human IL2.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

31.2 pg/mL-1000 pg/mL

Sensitivity

14 pg/mL