



## User's Manual

# CX3CL1 Human ELISA Kit



DEIA1401



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human CX3CL1 concentrations in Serum, Plasma, Cell Culture Supernatant and Urine.

### Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

0.25 ng/mL-60 ng/mL

### Sensitivity

300 pg/mL

