



User's Manual

CD105 Human ELISA Kit



DEIA1418



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human CD105 concentrations in serum, plasma, cell culture supernatants and urine.

Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

12.35 pg/mL-3000 pg/mL

Sensitivity

10 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com