



User's Manual

Rat Kit Ligand, Kitlg ELISA Kit

REF

DEIA142



5 x 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human DKK1 ELISA kit is for the quantitative determination of Human DKK1.

This ELISA kit contains the basic components required for the development of sandwich ELISAs. Each kit contains sufficient materials to run ELISAs on five 96-well plates.

Storage

Store at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

31.25 pg/mL-2,000 pg/mL

Detection Limit

31.25pg/mL

Sensitivity

31.25 pg/mL