



User's Manual

TRAIL Human ELISA Kit



DEIA1522



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human soluble TRAIL/Apo2L ELISA is to be used for the quantitative determination of soluble TNF-related apoptosis-inducing ligand (TRAIL) in Human serum, plasma, buffered solutions or cell culture medium.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

93.75 pg/mL-3000 pg/mL

Sensitivity

64 pg/mL