



User's Manual

Cytochrome-C Oxidase (Complex IV) Human Profiling ELISA Kit

REF

DEIA1523L



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Cytochrome-C oxidase (Complex IV) Human Profiling kit is an in vitro enzyme-linked immunosorbent assay (ELISA) for the comparison of cytochrome-C oxidase levels or profile in cell and tissue lysates.

Storage

Store all components at 4°C.

Precision

Intra-assay (CV%): 5.6%

Inter-assay (CV%): 10.9%

Detection Range

5 µg/mL-1000 µg/mL

Sensitivity

20 µg/ml

