



**User's Manual**

# **IL31 Human ELISA Kit (with color giving dyes)**



**DEIA1528**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human IL31 concentrations in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

7.8 pg/mL-500 pg/mL

### Sensitivity

1 pg/mL