



## User's Manual

# IL4 Mouse ELISA Kit



DEIA1530



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Mouse IL4 concentrations in cell and tissue lysates.

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

0.819 pg/mL-200 pg/mL

### Sensitivity

1 pg/mL