



## User's Manual

# IL6 Rat ELISA Kit



DEIA1540



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Rat IL6 concentrations in serum, plasma, cell culture supernatant and urine. This product is for research use only and is not intended for diagnostic use.

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

40.96 pg/mL-10000 pg/mL

### Sensitivity

30 pg/mL