



User's Manual

IL8 Pig ELISA Kit



DEIA1549



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Pig IL8 concentrations in serum, plasma and cell culture supernatants.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

8.19 pg/mL-2000 pg/mL

Sensitivity

10 pg/mL

