



User's Manual

Interferon gamma Human ELISA Kit



DEIA1551



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human Interferon gamma ELISA is to be used for the quantitative determination of interferon gamma in Human serum, plasma, buffered solutions or cell culture medium. The assay will recognize both natural and recombinant Human Interferon gamma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

12.5 pg/mL-400 pg/mL

Sensitivity

5 pg/mL