



User's Manual

Interferon gamma Human ELISA Kit



DEIA1552



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human Interferon gamma concentrations in serum, plasma, cell culture supernatant and urine.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

20.6 pg/mL-15000 pg/mL

Sensitivity

15 pg/mL