



User's Manual

Leptin Human ELISA Kit

REF

DEIA1560



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human Leptin concentrations in Serum, Plasma, Cell Culture Supernatant and Urine.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

2 pg/ml - 400 pg/ml

Sensitivity

2 pg/ml

