



**User's Manual**

# **Macrophage Inflammatory Protein 1 alpha / CCL3 Human ELISA Kit**

**REF**

**DEIA1566**



**96T**

**RUO**

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human Macrophage Inflammatory Protein 1 alpha/CCL3 concentrations in serum, plasma, cell culture supernatant and urine.

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

4.12 pg/mL-3000 pg/mL

### Sensitivity

6 pg/mL