



User's Manual

Macrophage Inflammatory Protein 1 alpha / CCL3 Mouse ELISA Kit

REF

DEIA1567



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Mouse Macrophage Inflammatory Protein 1 alpha/CCL3 concentrations in serum, plasma and cell culture supernatant.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

4.1 pg/mL-1000 pg/mL

Sensitivity

8 pg/mL

