



User's Manual

Osteoprotegerin Human ELISA Kit



DEIA1586



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human Osteoprotegerin concentrations in serum, plasma and cell culture supernatant and urine.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

1.23 pg/mL-900 pg/mL

Sensitivity

1 pg/mL

