



User's Manual

PAI1 + tPA Human ELISA



DEIA1588



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human Plasminogen Activator Inhibitor 1 (PAI-1) and Tissue plasminogen activator (tPA) concentrations in cell culture supernatants, tissue, plasma and milk.

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.

Detection Range

0.25 ng/mL - 16 ng/mL

Sensitivity

200 pg/mL

Specificity

Species: % Cross Reactivity

Canine: None

Bovine: None

Monkey: None

Mouse: None

Rat: None

Swine: None

Rabbit: 15