



User's Manual

Resistin Human ELISA Kit (with color giving dyes)



DEIA1599



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Human Resistin concentrations in cell culture supernatant, serum and plasma (EDTA, citrate and heparin).

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.

Detection Range

31 pg/mL-2000 pg/mL

Sensitivity

3.1 pg/mL