



**User's Manual**

# Retinol binding Protein 4 Human ELISA Kit



DEIA1603



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human Retinol binding Protein 4 concentrations in plasma, serum, saliva, urine and cell culture supernatants.

### Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

6.25 ng/mL-400 ng/mL

### Sensitivity

6 ng/mL