



User's Manual

Thyroxine (total T4) Human ELISA Kit



DEIA1617



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative measurement of Thyroxine (total T4) in Human serum and plasma.

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.

Sensitivity

0.004 µg/mL