



## User's Manual

# TIMP1 Human ELISA Kit



DEIA1618



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human TIMP1 concentrations in Serum, Plasma, Cell Culture Supernatant and Urine.

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

24.69 pg/mL-18000 pg/mL

### Sensitivity

40 pg/mL

